

Product Datasheet

PSPH Antibody (orb1246408)

Description	PSPH Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246408 (0.05 ug/ml) staining of Mous...
Immunogen	The immunogen for this antibody is: C-RQQVKDNAKWIYTD	
Target	PSPH	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 26kDa band observed in Human Brain (Hippocampus) and in Mouse brain lysates (calculated MW of 25.0kDa according to Human NP_004568.2 and 25.1kDa according to Mouse NP_598661.1). Recommended concentration: 0.05-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P78330	
NCBI	NP_004568.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 26kDa band observed in Human Brain (Hippocampus) and in Mouse brain lysates (calculated MW of 25.0kDa according to Human NP_004568.2 and 25.1kDa according to Mouse NP_598661.1). Recommended concentration: 0.05-1ug/ml.	
Expiration Date	12 months from date of receipt.	