

Product Datasheet

TMPRSS4 Antibody (orb1246484)

Description	TMPRSS4 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246484 (1 ug/ml) staining of Human B...
Immunogen	The immunogen for this antibody is: DSTRCNADDAYQGE	
Target	TMPRSS4	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 48kDa band observed in Human Bone Marrow lysates (calculated MW of 48.2kDa according to NP_063947.1 and of 47.6kDa according to NP_001077416.1). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9NRS4	
NCBI	NP_001077416.1 , NP_063947.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx 48kDa band observed in Human Bone Marrow lysates (calculated MW of 48.2kDa according to NP_063947.1 and of 47.6kDa according to NP_001077416.1). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	