

---

## Product Datasheet

### tfap2a Antibody (orb1246488)

Description	tfap2a Antibody
<b>Species/Host</b>	Goat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	The immunogen for this antibody is: C-NPNSHSETGSKAGDK
<b>Target</b>	tfap2a
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	500 ug/mL
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q8UVE6</a>
<b>NCBI</b>	<a href="#">NP_789829.1</a>
<b>Dilution Range</b>	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
<b>Expiration Date</b>	12 months from date of receipt.