
Product Datasheet

PRUNE2 Antibody (orb1246505)

Description

PRUNE2 Antibody


Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 orb1246505
(1 ug/ml)
staining of
Human B...

Immunogen

The immunogen for this antibody is: C-MESEKISEKQEE

Target

PRUNE2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 30kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 29.5kDa according to NP_620173.1). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

Uniprot ID
[Q8WUY3](#)
NCBI
[NP_620173.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 30kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 29.5kDa according to NP_620173.1). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.