



Product Datasheet

HSD17B3 antibody (PE) (orb124664)


Description	Rabbit polyclonal to HSD17B3 (PE).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human HSD17B3
Target	HSD17B3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124664 is the PE conjugated form of orb5476
Isotype	IgG
Clonality	Polyclonal
MW	34kDa
Uniprot ID	P37058
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of mouse brain lys...



Immunohistochemical analysis of formalin...