



Product Datasheet

CLCA1 Antibody (orb1246651)

Description	CLCA1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246651 (0.3 ug/ml) staining of Human...
Immunogen	The immunogen for this antibody is: PETPSPDETSAPC	
Target	CLCA1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1246651 (3.75 ug/ml) staining of para...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx. 100kDa band observed in Human Kidney and Lung lysates (calculated MW of 100kDa according to NP_001276.2). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	A8K7I4	
NCBI	NP_001276.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx. 100kDa band observed in Human Kidney and Lung lysates (calculated MW of 100kDa according to NP_001276.2). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	