

---

## Product Datasheet

### MTHFD2L Antibody (orb1246693)

**Description**

MTHFD2L Antibody


**Species/Host**

Goat

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

 orb1246693  
(0.1 ug/ml)  
staining of  
Human...

**Immunogen**

The immunogen for this antibody is: KHIQKEIQRGVES-C

**Target**

MTHFD2L

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Form/Appearance**

Liquid

**Concentration**

500 ug/mL

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 28kDa band observed in Human Testis lysates (calculated MW of 31.3kDa according to NP\_001004346.2). Recommended concentration: 0.1-0.3ug/ml.

**Clonality**

Polyclonal

**Uniprot ID**
[Q9H903](#)
**NCBI**
[NP\\_001004346.2](#)
**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 28kDa band observed in Human Testis lysates (calculated MW of 31.3kDa according to NP\_001004346.2). Recommended concentration: 0.1-0.3ug/ml.

**Expiration Date**

12 months from date of receipt.