

## Product Datasheet

# IFNAR1 (phospho-Ser535/Ser539) antibody (PE) (orb124672)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | IFNAR1 (phospho-Ser535/Ser539) antibody (PE)   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Rat   |
| <b>Conjugation</b>         | PE   |
| <b>Tested Applications</b> | FC, IF   |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human IFNAR1 around the phosphorylation site of Ser535+Ser539 |
| <b>Target</b>              | phospho-IFNAR1 (Ser535+Ser539)   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                     |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P17181</b>  |
| <b>Dilution Range</b>      | Flow-Cyt=2ug/tets,IF=1:100-500   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |