

Product Datasheet

IL17A antibody (PE) (orb124680)

Description	Rabbit polyclonal to IL17A (PE).
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from mouse IL-17 (85-150aa)
Target	IL17A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124680 is the PE conjugated form of orb10898
Isotype	IgG
Clonality	Polyclonal
MW	15kDa
Uniprot ID	Q62386
Expiration Date	12 months from date of receipt.