

Product Datasheet

TSPYL2 Antibody (orb1246830)

Description	TSPYL2 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246830 (0.3 ug/ml) staining of Jurka...
Immunogen	The immunogen for this antibody is: C-SDIDETIHDIKISD	
Target	TSPYL2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 79.4kDa according to NP_071400.1). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
MW	Approx 75 kDa	
Uniprot ID	Q9H2G4	
NCBI	NP_071400.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 79.4kDa according to NP_071400.1). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	