

Product Datasheet

TACR1 Antibody (orb1246881)

Description	TACR1 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: C-EIPGDSSDRYHEQ
Target	TACR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 36kDa band observed in Human Testis lysates (calculated MW of 35.7kDa according to NP_056542.1). Recommended concentration: 1-3ug/ml.
Clonality	Polyclonal
Uniprot ID	P25103
NCBI	NP_056542.1 , NP_001049.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 36kDa band observed in Human Testis lysates (calculated MW of 35.7kDa according to NP_056542.1). Recommended concentration: 1-3ug/ml.
Expiration Date	12 months from date of receipt.



250kDa
150kDa
100kDa



50kDa
37kDa
25kDa
20kDa
15kDa

orb1246881
(1 ug/ml)
staining of
Human T...