
Product Datasheet

SYNGR2 Antibody (orb1246889)

Description

SYNGR2 Antibody


Species/Host

Goat

Reactivity


Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


orb1246889
(1 ug/ml)
staining of
Human M...

Immunogen

The immunogen for this antibody is: CIYGEGYSNAHESKQ

Target

SYNGR2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 25-26kDa band observed in Human Muscle lysates (calculated MW of 24.8kDa according to NP_004701.1). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour.

Clonality

Polyclonal

Uniprot ID
[O43760](#)
NCBI
[NP_001307452.1](#), [NP_004701.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 25-26kDa band observed in Human Muscle lysates (calculated MW of 24.8kDa according to NP_004701.1). Recommended concentration: 0.3-1ug/ml. Primary incubation was 1 hour.

Expiration Date

12 months from date of receipt.