
Product Datasheet

SLC50A1 Antibody (orb1246977)

Description

SLC50A1 Antibody

Species/Host

Goat

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

The immunogen for this antibody is: C-EQDRNYWLLQT

Target

SLC50A1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments gave no signal but low background in Human Liver lysate at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Clonality

Polyclonal

Uniprot ID[Q9BRV3](#)**NCBI**[NP_001116309.1](#), [NP_061333.2](#), [NP_001116311.1](#)**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments gave no signal but low background in Human Liver lysate at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Expiration Date

12 months from date of receipt.