


Product Datasheet

TCF19 Antibody (orb1246980)

Description	TCF19 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: C-RSTAKAPSDTPAHE
Target	TCF19
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35 kDa band observed in 293 lysate (predicted MW of 41kDa according to NP_009040). Recommended for use at 1-2ug/ml.
Clonality	Polyclonal
Uniprot ID	Q9Y242
NCBI	NP_009040.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 35 kDa band observed in 293 lysate (predicted MW of 41kDa according to NP_009040). Recommended for use at 1-2ug/ml.
Expiration Date	12 months from date of receipt.



orb1246980 staining (1.5 ug/ml) of 293 l...