




Product Datasheet

SART1 Antibody (orb1246983)

Description	SART1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1246983 (0.03 ug/ml) staining of Hela...
Immunogen	The immunogen for this antibody is: GSSKKHRGEKEAA-C	
Target	SART1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 125kDa band observed in HeLa lysates (consistent with Vertegaal et al, below, calculated MW of 90.3 kDa according to NP_005137). Recommended concentration: 0.03-0.1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	O43290	
NCBI	NP_005137.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 125kDa band observed in HeLa lysates (consistent with Vertegaal et al, below, calculated MW of 90.3 kDa according to NP_005137). Recommended concentration: 0.03-0.1ug/ml.	
Expiration Date	12 months from date of receipt.	