



Product Datasheet

POLL Antibody (orb1247085)

Description	POLL Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-LGLPYREPAERDW	orb1247085 (5 ug/ml) staining of paraffi...
Target	POLL	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247085 (0.5 ug/ml) staining of A549 ...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 60kDa band observed in cell lysates of A549 and approx 60-70kDa in cell lysates of HepG2 (calculated MW of 63.5kDa according to Human NP_037406.1). Recommended concentration: 0.5-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9UGP5	
NCBI	NP_037406.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 60kDa band observed in cell lysates of A549 and approx 60-70kDa in cell lysates of HepG2 (calculated MW of 63.5kDa according to Human NP_037406.1). Recommended concentration: 0.5-1ug/ml.	
Expiration Date	12 months from date of receipt.	