
Product Datasheet

GOPC Antibody (orb1247096)

Description	GOPC Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The immunogen for this antibody is: C-LDDLHTLYHKKSY
Target	GOPC
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx. 55kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 50.5kDa according to NP_065132.1). Recommended concentration: 0.05-0.2ug/ml.Immunofluorescence: This product has been successfully used on HeLa cells.
Clonality	Polyclonal
Uniprot ID	Q9HD26
NCBI	NP_065132.1 , NP_001017408.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000.Western Blot:Approx. 55kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 50.5kDa according to NP_065132.1). Recommended concentration: 0.05-0.2ug/ml.Immunofluorescence: This product has been successfully used on HeLa cells.
Expiration Date	12 months from date of receipt.



orb1247096
(0.05 ug/ml)
staining of
Huma...