



Product Datasheet

PIK3C2A Antibody (orb1247103)

Description	PIK3C2A Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	orb1247103 (2 ug/ml) staining of paraffi...
Immunogen	The immunogen for this antibody is: C-TVKQYQLTAATYL	
Target	PIK3C2A	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Testis and Molt-4 lysates at up to 1ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows cytoplasm staining of cytotrophoblasts and mesenchymal cells. Recommended concentration: 2-4ug/ml. Paraffin embedded Human Heart. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	orb1247103 (5 ug/ml) staining of paraffi...
Uniprot ID	O00443	
NCBI	NP_002636.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Testis and Molt-4 lysates at up to 1ug/ml. Immunohistochemistry: In paraffin embedded Human Placenta shows cytoplasm staining of cytotrophoblasts and mesenchymal cells. Recommended concentration: 2-4ug/ml. Paraffin embedded Human Heart. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	