
Product Datasheet

PGAM1/PGAM2/PGAM4 Antibody (orb1247114)

Description

PGAM1/PGAM2/PGAM4 Antibody

Species/Host

Goat

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

The immunogen for this antibody is: C-KAMEAVAAQGKAKK

Target

PGAM1 / PGAM2 / PGAM4

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 28kDa band observed in Human, Mouse, Rat and Pig Liver lysates (calculated MW of 28.8kDa according to NP_002620.1; NP_000281.2 and of NP_001025062.1). The additional 30kDa band observed in the rodent lysates likely represent the acetylated form. Recommended concentration: 0.05-0.15ug/ml.

Clonality


Polyclonal

Uniprot ID
[P18669](#)
NCBI
[NP_000281.2](#), [NP_001025062.1](#), [NP_002620.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 28kDa band observed in Human, Mouse, Rat and Pig Liver lysates (calculated MW of 28.8kDa according to NP_002620.1; NP_000281.2 and of NP_001025062.1). The additional 30kDa band observed in the rodent lysates likely represent the acetylated form. Recommended concentration: 0.05-0.15ug/ml.

Expiration Date

12 months from date of receipt.



orb1247114
(0.05 ug/ml)
staining of
Huma...