





Product Datasheet

PCBP4 Antibody (orb1247137)

Description	PCBP4 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1247137 (0.03 ug/ml) staining of HeLa...
Immunogen	The immunogen for this antibody is: C-NGSKKAERQKFSKY	
Target	PCBP4	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	 HEK293 overexpressing PCBP4 and probed w...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx. 40kDa band observed in lysates of cell line HeLa (calculated MW of 41.4kDa according to NP_127502.1 and NP_127503.1). In transfected HEK293 transiently expressing PCBP4 a band of approx. 45kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	P57723	
NCBI	NP_127503.1 , NP_127501.1 , NP_127502.1 , NP_065151.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx. 40kDa band observed in lysates of cell line HeLa (calculated MW of 41.4kDa according to NP_127502.1 and NP_127503.1). In transfected HEK293 transiently expressing PCBP4 a band of approx. 45kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	