

Product Datasheet

ARHGEF18 Antibody (orb1247154)

Description	ARHGEF18 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1247154 (0.5 ug/ml) staining of lysat...
Immunogen	The immunogen for this antibody is: C-SAKEDASKEDVIFF	
Target	ARHGEF18	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 140+90kDa band observed in lysates of cell line Daudi and Caco-2 (calculated MW of 114kDa according to NP_056133.2 and 131kDa according to NP_001124427.1). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q6ZSZ5	
NCBI	NP_001124427.1 , NP_056133.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 140+90kDa band observed in lysates of cell line Daudi and Caco-2 (calculated MW of 114kDa according to NP_056133.2 and 131kDa according to NP_001124427.1). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	