





---

## Product Datasheet

### ABCC5 Antibody (orb1247251)

<b>Description</b>	ABCC5 Antibody	
<b>Species/Host</b>	Goat	
<b>Reactivity</b>	Human	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	ELISA, IHC, WB	
<b>Immunogen</b>	The immunogen for this antibody is: KDIDIGKEYIIP-C	orb1247251 (0.5 ug/ml) staining of Human...
<b>Target</b>	ABCC5	
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	500 ug/mL	
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Note</b>	For research use only	orb1247251 (4 ug/ml) staining of paraffi...
<b>Application notes</b>	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 150kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 161kDa according to NP_005679.2). Recommended concentration: 0.5-1.5ug/ml. Immunohistochemistry: In paraffin embedded Human Kidney shows membranous staining in DCT. Paraffin embedded Human Liver. Recommended concentration: 4-6ug/ml.	
<b>Clonality</b>	Polyclonal	orb1247251 (5 ug/ml) staining of paraffi...
<b>Uniprot ID</b>	<a href="#">O15440</a>	
<b>NCBI</b>	<a href="#">NP_001018881.1</a> , <a href="#">NP_005679.2</a>	
<b>Dilution Range</b>	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx. 150kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 161kDa according to NP_005679.2). Recommended concentration: 0.5-1.5ug/ml. Immunohistochemistry: In paraffin embedded Human Kidney shows membranous staining in DCT. Paraffin embedded Human Liver. Recommended concentration: 4-6ug/ml.	
<b>Expiration Date</b>	12 months from date of receipt.	