

Product Datasheet

IFT57 Antibody (orb1247256)

Description	IFT57 Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-EQPQEYDDPNATISN	orb1247256 (1 ug/ml) staining of Mouse B...
Target	IFT57	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 49kDa band observed in Human and Mouse Brain lysates (calculated MW of 49.1kDa according to Human NP_060480.1 and of 48.8kDa according to Mouse NP_082956.2). Recommended concentration: 1-3ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9NWB7	
NCBI	NP_060480.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 49kDa band observed in Human and Mouse Brain lysates (calculated MW of 49.1kDa according to Human NP_060480.1 and of 48.8kDa according to Mouse NP_082956.2). Recommended concentration: 1-3ug/ml.	
Expiration Date	12 months from date of receipt.	