





Product Datasheet

MAOA Antibody (orb1247259)

Description	MAOA Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	orb1247259 (0.3 ug/ml) staining of Human...
Immunogen	The immunogen for this antibody is: C-DAPWEAQHADKWDK	
Target	MAOA	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247259 (5 ug/ml) staining of paraffi...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 70kDa band observed in Human Heart lysates (calculated MW of 60.0kDa according to NP_000231.1). The user is warned for non-specific background in Western blot. Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Placenta and Small Intestine. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	orb1247259 (5 ug/ml) staining of paraffi...
Uniprot ID	P21397	
NCBI	NP_000231.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 70kDa band observed in Human Heart lysates (calculated MW of 60.0kDa according to NP_000231.1). The user is warned for non-specific background in Western blot. Recommended concentration: 0.1-0.3ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Placenta and Small Intestine. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	