

Product Datasheet

UBE2V2 Antibody (orb1247261)

Description	UBE2V2 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: AVSTGVKVP RN F-C
Target	UBE2V2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 20kDa band seen in human brain lysate. Recommended for use at 0.2-1ug/ml.
Clonality	Polyclonal
Uniprot ID	Q15819
NCBI	NP_003341.1 , NP_954673.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 20kDa band seen in human brain lysate. Recommended for use at 0.2-1ug/ml.
Expiration Date	12 months from date of receipt.



orb1247261 staining (0.5 ug/ml) of human...