
Product Datasheet

FBXL12 Antibody (orb1247523)

Description

FBXL12 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Immunogen

The immunogen for this antibody is: RACPKESMDWWM

Target

FBXL12

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Human Lung shows staining of ciliated epithelial cells in the bronchioles. Recommended concentration, 4-6ug/ml.

Clonality


Polyclonal

Uniprot ID
[Q9NXX8](#)
NCBI
[NP_060173.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1ug/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Human Lung shows staining of ciliated epithelial cells in the bronchioles. Recommended concentration, 4-6ug/ml.

Expiration Date

12 months from date of receipt.



orb1247523 (4 ug/ml) staining of paraffi...