


Product Datasheet

FAF2 Antibody (orb1247543)

Description	FAF2 Antibody	
Species/Host	Goat	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1247543 staining (0.3 ug/ml) of Jurka...
Immunogen	The immunogen for this antibody is: C-HTEVLVQDLTDE	
Target	FAF2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247543 (1 ug/ml) staining of Human (...)
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in lysates of Human and Mouse Thymus and of cell line Jurkat (calculated MW of 52.6kDa according to NP_055428.1). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q96CS3	
NCBI	NP_055428.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 55kDa band observed in lysates of Human and Mouse Thymus and of cell line Jurkat (calculated MW of 52.6kDa according to NP_055428.1). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	