


Product Datasheet

TPD52L1 Antibody (orb1247627)

Description	TPD52L1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	The immunogen for this antibody is: C-EPLQGTDEDAVASAD	
Target	TPD52L1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	orb1247627 (0.1 ug/ml) staining of Human...
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 26kDa band observed in Human Breast cancer lysates (calculated MW of 22.5kDa according to NP_003278.1). Recommended concentration: 0.1-0.3ug/ml. See also our D52 and D54 products.	
Clonality	Polyclonal	
Uniprot ID	Q16890	
NCBI	NP_001003397.1 , NP_001003395.1 , NP_001003396.1 , NP_003278.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 26kDa band observed in Human Breast cancer lysates (calculated MW of 22.5kDa according to NP_003278.1). Recommended concentration: 0.1-0.3ug/ml. See also our D52 and D54 products.	
Expiration Date	12 months from date of receipt.	