

---

## Product Datasheet

### ACTL6A/ACTL6B Antibody (orb1247745)

**Description**

ACTL6A/ACTL6B Antibody


 250kDa  
150kDa

**Species/Host**

Goat

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

 orb1247745  
staining (1  
ug/ml) of  
HeLa ly...

**Immunogen**

The immunogen for this antibody is: YEEGGKQCVERKCP

**Target**

ACTL6A / ACTL6B

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Form/Appearance**

Liquid

**Concentration**

500 ug/mL

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 48kDa band observed in HeLa lysates (predicted MW of 49kDa according to NP\_057272.1, NP\_004292.1). Recommended for use at 1-3ug/ml. Primary incubation was for 1 hour.

**Clonality**

Polyclonal

**Uniprot ID**
[O94805](#)
**NCBI**
[NP\\_057272.1](#), [NP\\_817126.1](#), [NP\\_004292.1](#)
**Dilution Range**

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 48kDa band observed in HeLa lysates (predicted MW of 49kDa according to NP\_057272.1, NP\_004292.1). Recommended for use at 1-3ug/ml. Primary incubation was for 1 hour.

**Expiration Date**

12 months from date of receipt.