

Product Datasheet

ARL4A Antibody (orb1247770)

Description	ARL4A Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: GNGLSDQTSILSN-C
Target	ARL4A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:1000.Western Blot:Approx 23kDa band seen in Human Testis lysate. Recommended for use at 1-3ug/ml.
Clonality	Polyclonal
Uniprot ID	P40617
NCBI	P_997625.1 , NP_005729.1 , NP_001032241.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000.Western Blot:Approx 23kDa band seen in Human Testis lysate. Recommended for use at 1-3ug/ml.
Expiration Date	12 months from date of receipt.



orb1247770 staining (2 ug/ml) of Human T...