





---

## Product Datasheet

### APH1A Antibody (orb1247790)

<b>Description</b>	APH1A Antibody	
<b>Species/Host</b>	Goat	
<b>Reactivity</b>	Human	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	ELISA, IHC, WB	orb1247790 (0.05 ug/ml) staining of Huma...
<b>Immunogen</b>	The immunogen for this antibody is: C-HVTDRSDARLQYG	
<b>Target</b>	APH1A	
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	500 ug/mL	
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Note</b>	For research use only	orb1247790 (5 ug/ml) staining of paraffi...
<b>Application notes</b>	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 26.8kDa according to NP_057106.2). Recommended concentration: 0.05-0.1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Placenta and Kidney. Recommended concentration: 5ug/ml.	
<b>Clonality</b>	Polyclonal	
<b>Uniprot ID</b>	<a href="#">Q96BI3</a>	orb1247790 (5 ug/ml) staining of paraffi...
<b>NCBI</b>	<a href="#">NP_001071096.1</a> , <a href="#">NP_057106.2</a>	
<b>Dilution Range</b>	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 26kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 26.8kDa according to NP_057106.2). Recommended concentration: 0.05-0.1ug/ml. Primary incubation was 1 hour. Immunohistochemistry: Paraffin embedded Human Placenta and Kidney. Recommended concentration: 5ug/ml.	
<b>Expiration Date</b>	12 months from date of receipt.	