

Product Datasheet

SNTA1 Antibody (orb1247889)

Description

SNTA1 Antibody



Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1247889
(0.01 ug/ml)
staining of
Huma...

Immunogen

The immunogen for this antibody is: ASGRRAPRTGLLE-C

Target

SNTA1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 60kDa band observed in Human and Mouse Muscle lysates and in rodent Heart lysates (calculated MW of 53.9kDa according to Human NP_003089.1). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q13424](#)

NCBI

[NP_003089.1](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 60kDa band observed in Human and Mouse Muscle lysates and in rodent Heart lysates (calculated MW of 53.9kDa according to Human NP_003089.1). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.