



Product Datasheet

SH3GLB1 Antibody (orb1247894)

Description	SH3GLB1 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	The immunogen for this antibody is: C-QKGKVPITYLELLN	orb1247894 staining (0.3 ug/ml) of Jurka...
Target	SH3GLB1	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247894 (2.5 ug/ml) staining of paraf...
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 38kDa band observed in lysates of Human T cell line Jurkat and in Human Lung and Placenta (calculated MW of 40.8kDa according to NP_057093.1). Recommended concentration: 0.3-1.0 ug/ml. A minor band of unknown identity was also consistently observed at 50kDa in Jurkat and Placenta. This band was successfully blocked by incubation with the immunising peptide. Immunohistochemistry: Paraffin embedded Human Colon. Recommended concentration: 2.5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9Y371	
NCBI	NP_057093.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 38kDa band observed in lysates of Human T cell line Jurkat and in Human Lung and Placenta (calculated MW of 40.8kDa according to NP_057093.1). Recommended concentration: 0.3-1.0 ug/ml. A minor band of unknown identity was also consistently observed at 50kDa in Jurkat and Placenta. This band was successfully blocked by incubation with the immunising peptide. Immunohistochemistry: Paraffin embedded Human Colon. Recommended concentration: 2.5ug/ml.	
Expiration Date	12 months from date of receipt.	