


---

## Product Datasheet

### SEC61A1 Antibody (orb1247896)

<b>Description</b>	SEC61A1 Antibody	
<b>Species/Host</b>	Goat	
<b>Reactivity</b>	Human	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	ELISA, WB	orb1247896 (0.3 ug/ml) staining of Human...
<b>Immunogen</b>	The immunogen for this antibody is: C-KEQSEVGSMSGALLF	
<b>Target</b>	SEC61A1	
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	500 ug/mL	
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Note</b>	For research use only	
<b>Application notes</b>	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 48kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 52.3kDa according to NP_037468.1). Recommended concentration: 0.3-1ug/ml.	orb1247896 (3.75 ug/ml) staining of para...
<b>Clonality</b>	Polyclonal	
<b>Uniprot ID</b>	<a href="#">P61619</a>	
<b>NCBI</b>	<a href="#">NP_037468.1</a>	
<b>Dilution Range</b>	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Approx 48kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 52.3kDa according to NP_037468.1). Recommended concentration: 0.3-1ug/ml.	
<b>Expiration Date</b>	12 months from date of receipt.	