



Product Datasheet

FTCD Antibody (orb1247976)

Description	FTCD Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1247976 (0.05 ug/ml) staining of Huma...
Immunogen	The immunogen for this antibody is: SQLVECVPNFSEGKNQ	
Target	FTCD	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1247976 (3.8 ug/ml) staining of paraf...
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 60kDa band observed in Human Liver lysates (calculated MW of 58.9kDa according to NP_006648.1. Recommended concentration: 0.05-0.2ug/ml. Primary incubation was 1 hour.	
Clonality	Polyclonal	
Uniprot ID	O95954	
NCBI	NP_006648.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 60kDa band observed in Human Liver lysates (calculated MW of 58.9kDa according to NP_006648.1. Recommended concentration: 0.05-0.2ug/ml. Primary incubation was 1 hour.	
Expiration Date	12 months from date of receipt.	