

Product Datasheet

SP140 Antibody (orb1248130)

Description	SP140 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: EGQNLQEQVCPEP
Target	SP140
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 140kDa band observed in Human Spleen lysates (calculated MW of 98.2kDa according to NP_009168.4). Recommended concentration: 1-3ug/ml.
Clonality	Polyclonal
MW	Approx 140kDa
Uniprot ID	Q13342
NCBI	NP_001265382.1 , NP_001265380.1 , NP_001005176.1 , NP_009168.4 , NP_001265381.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 140kDa band observed in Human Spleen lysates (calculated MW of 98.2kDa according to NP_009168.4). Recommended concentration: 1-3ug/ml.
Expiration Date	12 months from date of receipt.



orb1248130
(2 ug/ml)
staining of
Human S...