


Product Datasheet

Acacb Antibody (orb1248269)

Description	Acacb Antibody	 250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa 15kDa
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	orb1248269 (5 ug/ml) staining of paraffi...
Immunogen	The immunogen for this antibody is: C-TLNVQESGRVKYIKR	
Target	Acacb	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 277kDa according to NP_001084.3. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. Immunohistochemistry: Paraffin embedded Human Adrenal Gland. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	E9Q4Z2	
NCBI	NP_598665.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. Western Blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 277kDa according to NP_001084.3. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. Immunohistochemistry: Paraffin embedded Human Adrenal Gland. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	