

Product Datasheet

ANAPC11 Antibody (orb1248304)

Description	ANAPC11 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	orb1248304 (1 ug/ml) staining of Cerebra...
Immunogen	The immunogen for this antibody is: PALDQETSSLLRCTS	
Target	ANAPC11	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 19kDa band observed in Human Brain (Cerebral and Frontal Cortex, Hippocampus) lysates (calculated MW of 20.6kDa according to NP_001002244.1). Recommended concentration: 0.3-1ug/ml.	
Clonality	Polyclonal	
MW	Approx 19kDa	
Uniprot ID	Q9NYG5	
NCBI	NP_001002244.1 , NP_057560.8	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 19kDa band observed in Human Brain (Cerebral and Frontal Cortex, Hippocampus) lysates (calculated MW of 20.6kDa according to NP_001002244.1). Recommended concentration: 0.3-1ug/ml.	
Expiration Date	12 months from date of receipt.	