

Product Datasheet

Nalp12/Nlrp12 Antibody (orb1248314)

Description

Nalp12/Nlrp12 Antibody


Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 orb1248314
(1 ug/ml)
staining of
Mouse S...

Immunogen

 The immunogen for this antibody is: C-
KPGTPTFKVPANQR

Target

Nalp12 / Nlrp12

Preservatives

 Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium
azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot
and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

 Aliquot and store at -20°C. Minimize freezing and
thawing.

Note

For research use only

Application notes

 Peptide ELISA: antibody detection limit dilution
1:128000. Western Blot: Approx 100kDa band observed
in Mouse Liver and Spleen lysates (calculated MW of
119kDa according to NP_001028603.1). The observed
molecular weight corresponds to earlier findings in
literature with different antibodies (Jéru et al, Proc Natl
Acad Sci U S A. 2008 Feb 5;105(5):1614-9. PMID:
18230725). Recommended concentration: 1-3ug/ml.

Clonality

Polyclonal

MW

Approx 100kDa

Uniprot ID
[E9Q5R7](#)
NCBI
[NP_001028603.1](#)
Dilution Range

 Peptide ELISA: antibody detection limit dilution
1:128000. Western Blot: Approx 100kDa band observed
in Mouse Liver and Spleen lysates (calculated MW of
119kDa according to NP_001028603.1). The observed
molecular weight corresponds to earlier findings in
literature with different antibodies (Jéru et al, Proc Natl
Acad Sci U S A. 2008 Feb 5;105(5):1614-9. PMID:
18230725). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.