


Product Datasheet

PIM2 Antibody (orb1248403)

Description	PIM2 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	HEK293 lysate (10 ug protein in RIPA buf...
Immunogen	The immunogen for this antibody is: C-KDREAFAEYRLGP	
Target	PIM2	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	orb1248403 (5 ug/ml) staining of paraffi...
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: In transfected HEK293 transiently expressing Human PIM2 (myc and DYKDDDDK tagged), a band of approx. 38kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration: 0.5-1ug/ml. Immunohistochemistry: Paraffin embedded Human Heart. Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q9P1W9	
NCBI	NP_006866.2	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: In transfected HEK293 transiently expressing Human PIM2 (myc and DYKDDDDK tagged), a band of approx. 38kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration: 0.5-1ug/ml. Immunohistochemistry: Paraffin embedded Human Heart. Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	