
Product Datasheet

HWP1 Antibody (orb1248413)

Evolution Biocomponents.
Description

HWP1 Antibody

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The immunogen for this antibody is: C-QPSIPAGETSPAVPK
Target	HWP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product.
Clonality	Polyclonal
Uniprot ID	P46593
NCBI	XP_712998.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product.
Expiration Date	12 months from date of receipt.