
Product Datasheet

Contactin 4/Big-2 Antibody (orb1248545)

Description

Contactin 4/Big-2 Antibody



250kDa
150kDa
100kDa

Species/Host

Goat

Reactivity


Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


50kDa
37kDa
25kDa
20kDa
15kDa

orb1248545
(0.5 ug/ml)
staining of
Mouse...

Immunogen

The immunogen for this antibody is: C-RRADGKPIARK

Target

Contactin 4 / Big-2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 150kDa band observed in Mouse Brain lysates and not in the same lysates from the KO mice (calculated MW of 113kDa according to NP_001103219.1). Recommended concentration: 0.5-1.5ug/ml.

Clonality

Polyclonal

MW

Approx 150kDa

Uniprot ID
[Q81WV2](#)
NCBI
[NP_766592.2](#), [NP_001103219.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000. Western Blot: Approx 150kDa band observed in Mouse Brain lysates and not in the same lysates from the KO mice (calculated MW of 113kDa according to NP_001103219.1). Recommended concentration: 0.5-1.5ug/ml.

Expiration Date

12 months from date of receipt.