

Product Datasheet

Cyb5r4 Antibody (orb1248701)

Description	Cyb5r4 Antibody
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen for this antibody is: RCQLEKLALREKRFD
Target	Cyb5r4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Rat Spleen lysates (calculated MW of 55kDa according to NP_596918.2 and 58.8kDa according to Q68EJ0-1). Recommended concentration: 0.1-0.3ug/ml.
Clonality	Polyclonal
MW	Approx 65 kDa
Uniprot ID	Q3TDX8
NCBI	NP_077157.2
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 65kDa band observed in Rat Spleen lysates (calculated MW of 55kDa according to NP_596918.2 and 58.8kDa according to Q68EJ0-1). Recommended concentration: 0.1-0.3ug/ml.
Expiration Date	12 months from date of receipt.



orb1248701
(0.1 ug/ml)
staining of
Rat S...