
Product Datasheet

AGTPBP1 Antibody (orb1249018)

Description

AGTPBP1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

The immunogen for this antibody is: C-TSPLEYNLPS

Target

AGTPBP1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 150kDa band observed in lysates of cell line HeLa (calculated MW of 134kDa according to NP_056054.2). Recommended concentration: 1-3ug/ml.

Clonality


Polyclonal

Uniprot ID
[Q9UPW5](#)
NCBI
[NP_056054.2](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 150kDa band observed in lysates of cell line HeLa (calculated MW of 134kDa according to NP_056054.2). Recommended concentration: 1-3ug/ml.

Expiration Date

12 months from date of receipt.


 orb1249018
(1 ug/ml)
staining of
HeLa ly...