



Product Datasheet

BIRC3 Antibody (orb1249038)

Description	BIRC3 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen for this antibody is: C-SQPTFPSSVT
Target	BIRC3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 68.4kDa according to NP_001156.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil. Recommended concentration: 5ug/ml.
Clonality	Polyclonal
Uniprot ID	Q13489
NCBI	NP_001156.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 68.4kDa according to NP_001156.1). Recommended concentration: 0.1-0.3ug/ml. Immunohistochemistry: Paraffin embedded Human Tonsil. Recommended concentration: 5ug/ml.
Expiration Date	12 months from date of receipt.



orb1249038
(0.1 ug/ml)
staining of
HeLa ...



orb1249038 (5
ug/ml)
staining of
paraffi...