


Product Datasheet

TRIM71 Antibody (orb1249149)

| | | |
|----------------------------|--|---|
| Description | TRIM71 Antibody |  |
| Species/Host | Goat | |
| Reactivity | Human, Mouse | |
| Conjugation | Unconjugated | |
| Tested Applications | ELISA, IHC, WB | |
| Immunogen | The immunogen for this antibody is: C-RQVKAKSLYLQVE | orb1249149 (0.3 ug/ml) staining of HEK29... |
| Target | TRIM71 | |
| Preservatives | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing. | |
| Form/Appearance | Liquid | |
| Concentration | 500 ug/mL | |
| Storage | Aliquot and store at -20°C. Minimize freezing and thawing. | |
| Note | For research use only | |
| Application notes | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 100kDa band observed in lysates of the cell lines HEK293 and NIH3T3 (calculated MW of 93.4kDa according to Human NP_001034200.1 and 92.1kDa according to Mouse NP_001035968.1). The observed molecular weight corresponds to findings with EB09902, a product of different design reacting to a different part of the same protein. Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 3.75ug/ml. |  |
| Clonality | Polyclonal | orb1249149 (3.75 ug/ml) staining of para... |
| Uniprot ID | Q2Q1W2 | |
| NCBI | NP_001034200.1 | |
| Dilution Range | Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Approx 100kDa band observed in lysates of the cell lines HEK293 and NIH3T3 (calculated MW of 93.4kDa according to Human NP_001034200.1 and 92.1kDa according to Mouse NP_001035968.1). The observed molecular weight corresponds to findings with EB09902, a product of different design reacting to a different part of the same protein. Recommended concentration: 0.3-1ug/ml. Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 3.75ug/ml. | |
| Expiration Date | 12 months from date of receipt. | |