

---

## Product Datasheet

### ACPP Antibody (orb1249269)

Description	ACPP Antibody	
<b>Species/Host</b>	Goat	
<b>Reactivity</b>	Human	orb1249269 (0.01 ug/ml) staining of Huma...
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	ELISA, IHC, WB	
<b>Immunogen</b>	The immunogen for this antibody is: C-NHMKRATQIPSYKK	
<b>Target</b>	ACPP	
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	500 ug/mL	
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.	orb1249269 (3.8 ug/ml) staining of paraf...
<b>Note</b>	For research use only	
<b>Application notes</b>	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 48.3kDa according to NP_001127666.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows reticulate cytoplasm staining, more intens towards the lumen of the gland. Recommended concentration, 3-5ug/ml.	
<b>Clonality</b>	Polyclonal	
<b>Uniprot ID</b>	<a href="#">P15309</a>	
<b>NCBI</b>	<a href="#">NP_001090.2</a> , <a href="#">NP_001127666.1</a>	
<b>Dilution Range</b>	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 48.3kDa according to NP_001127666.1). Recommended concentration: 0.01-0.03ug/ml. Immunohistochemistry: In paraffin embedded Human Prostate shows reticulate cytoplasm staining, more intens towards the lumen of the gland. Recommended concentration, 3-5ug/ml.	
<b>Expiration Date</b>	12 months from date of receipt.	