

Product Datasheet

SLC46A1 Antibody (orb1249303)

Description

SLC46A1 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb1249303
(0.3 ug/ml)
staining of
Human...

Immunogen

The immunogen for this antibody is: C-LKEPKSTRLFTFRH

Target

SLC46A1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 50kDa band observed in Human Kidney lysates (calculated MW of 49.8kDa according to NP_542400.2). Recommended concentration: 0.3-1ug/ml.

Clonality

Polyclonal

Uniprot ID

[Q96NT5](#)

NCBI

[NP_542400.2](#)

Dilution Range

Peptide ELISA: antibody detection limit dilution 1:64000. Western Blot: Approx 50kDa band observed in Human Kidney lysates (calculated MW of 49.8kDa according to NP_542400.2). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.



250kDa
150kDa
100kDa



50kDa
37kDa
25kDa
20kDa
15kDa