

Product Datasheet

BCL2L15 Antibody (orb1249339)

Description	BCL2L15 Antibody	
Species/Host	Goat	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC	orb1249339 (5 ug/ml) staining of paraffi...
Immunogen	The immunogen for this antibody is: C-SPTLQVASRN	
Target	BCL2L15	
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Form/Appearance	Liquid	
Concentration	500 ug/mL	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in lysates of cell lines A549, HEK293 and Jurkat gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 5ug/ml.	
Clonality	Polyclonal	
Uniprot ID	Q5TBC7	
NCBI	NP_001010922.1	
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. Western Blot: Preliminary experiments in lysates of cell lines A549, HEK293 and Jurkat gave no specific signal but low background (at antibody concentration up to 1ug/ml). Immunohistochemistry: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 5ug/ml.	
Expiration Date	12 months from date of receipt.	