

Product Datasheet

LPPR3 Antibody (orb1249344)

Description

LPPR3 Antibody

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The immunogen for this antibody is: C-TQRGHDSVYQQN
Target	LPPR3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Liquid
Concentration	500 ug/mL
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Clonality	Polyclonal
Uniprot ID	Q6T4P5
NCBI	NP_079164.1
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Expiration Date	12 months from date of receipt.